

**INSTRUCTIONS**  
**for**  
**MODEL NO. HT-300/500**  
**HARDNESS TESTER**  
**VERSION 1.0**

KEY INTERNATIONAL, INC.  
480 ROUTE 9  
ENGLISHTOWN, NJ 07726  
Phone: (732) 536-9700 Fax: (732) 972-2630

### Table of Contents

General Description .....	4
Introduction .....	4
Specifications .....	4
Method of Operation .....	5
Putting into Operation.....	5
Wiring.....	5
Inserting the Chart Paper .....	6
Operating Keyboard.....	7
Operating controls .....	7
Power Button.....	7
Restart Switch .....	7
Stop Button.....	7
Start Button .....	7
Printer .....	7
Adjustments .....	8
Limit Switch .....	8
Unit .....	8
Switch for units .....	8
Adjustment.....	8
Starting the test program .....	8
Zero balance .....	8
Calibration.....	9
Explanation of error information.....	9
CONNECTION POSSIBILITIES .....	10

Connection of an external start-key .....10

    Connection to arallel-interface (94-00442) (OPTION).....10

Plug connections.....10

    a) 5-pin socket: .....10

Spare Parts List.....16

## General Description


### Introduction

The PTB 300 is a high performance hardness tester, which can accurately measure the tensile strength of tablets, ampoules, food and other articles, either in single or in series operation. The tensile strength can be printed out in Kilopond (Kp), Newton (N) or Strong Cobb (SC) units, according to the operator's specifications.

### Specifications

Sample Diameter	4 - 40 mm
Jaw Stroke	4 - 40 mm (externally adjustable)
Chart Speed	3 mm/second
Force Increase	2 Kp/second
Measuring Range	HT-300: 0.5 - 30 Kp HT-500: 0.8 - 50 Kp , (convertable to Newtons (N), Kilo Ponds (Kp), or Strong Cobbs (Sc) )
Measured Value	Printed on the built in printer
Resolution	0.1 Kp
Accuracy	+/- 1% Full Scale
Measuring Jaws	Interchangable for different sample forms (tablets, ampoules, etc...)
Dimensions	(W) 530 mm x (D) 240 mm x (H) 180 mm
Power Supply	Switch Selectable, 110/220 volt, 50/60 Hz., Single Phase
Input	External Start Switch (optional)
Output	RS-232 or BCD Parallel

## Method of Operation



**Figure 1: Method of Operation**

After pressing the START-key, the zero-value of the load cell is determined and stored in the microprocessor. If the value is within the permitted range, the motor, controlled by the regulator-unit and the motor-amplifier, causes the plunger to move at a constant speed towards the specimen. As soon as the plunger touches the specimen and produces a force on the load cell, the microprocessor switches the regulator on so that a linearly increasing force is produced via the control circuit (load cell-regulator-motor-specimen) until the specimen breaks. The decreasing signal from the load cell indicates to the microprocessor, which subtracts from this value the zero-value (automatic zero-point correction) and prints out the calculated value in the desired units (N, kp, SC).

The spring-unit in front of the load cell provides a smooth transition from the constant velocity phase of the plunger to the linear power increase phase.


## Putting into Operation

**NOTE: Before putting the HT-300 into operation, check the voltage setting switch on the back of the unit (220 V resp. 115 V) and make sure that the operating voltage corresponds to the local power supply.**

## Wiring

- connect the signal cables to any additional units
- connect the signal cable from the external start-key (if fitted)

- connect power supply line.



**Figure 2: Back of Unit**

1. Output 25-pin, plug female type D (parallel-interface)
2. Power input with fuse (250 mA at 220 Volt, 500 mA at 115 Volt)
3. Voltage selector switch 220 V / 115 V
4. Waste Container: Under the container is the switch for units selection, the potentiometer for zero adjustment and the potentiometer for balancing the load cell
5. Under the Waste container and behind the access door is the switch for units selection, the potentiometer for zero adjustment and the potentiometer for balancing the load cell
6. Motor fuse 400 mA
7. Input socket plug (optional external start switch)

### **Inserting the Chart Paper**

The instrument should never be operated without paper.

To insert the paper open the cover (1). Slide in a roll of paper (2), with the dim side up, onto the spindle (3).

Then feed the paper into the printer (4) and while turning the paper transport wheel (5) lead the paper through the printer.

## **Operating Keyboard**

### **Operating controls**

#### **Power Button**

Switch on the instrument by pressing the POWER button (1). When pressed the switch should illuminate to indicate that the power has been switched on.

#### **Restart Switch**

Using the thumb wheel selector switch (4) the operator can set the machine to automatically start a series of tests. When set to position 0 the automatic testing option is turned off. Selecting positions from 1 to 9 will activate the automatic operation mode. In the automatic mode the operator places a sample for testing between the jaws and presses the start button. After the test is completed the jaws will move away and automatically start the next test after the number of seconds as set by the restart switch. If any errors occur, see errors section in this manual, the automatic operation will stop and can be restarted by pressing the start button. After all the samples have been tested press the stop button or set the restart switch to 0 and the automatic operation will be stopped.

#### **Stop Button**

By pressing the STOP button the testing procedure can be interrupted at any time.

#### **Start Button**

- a) By pressing the start button (6) the single measuring procedure can be started (RESTART switch (4) must be set to O). The START button is only operative when the movable plunger is in neutral position.
- b) With the START button a series of measurements can be started (RESTART switch (4) must not be set to O).

#### **Printer**

When a printout is terminated and the printing paper has stopped moving, or when the paper is torn, move the release-lever (2) in direction of the arrow and remove the paper from the printer. The release lever must not be touched during operation.

The printer (3) should never be operated without paper, otherwise a clean printout can not be guaranteed.

## **Adjustments**

### **Limit Switch**

By means of the limit switch the jaw stroke can be adjusted from approximately 4 - 40 mm.

### **Unit**

Using this switch following units can be selected as required: Kp, N, SC.

### **Switch for units**

Units Selection Switch

Potentiometer "GAIN"

Potentiometer "ZERO"

### **Adjustment**

A test program is incorporated in the PTB 300 which allows to check and if necessary to adjust the zero-point- and amplification adjustment without the need of special measuring equipment.

### **Starting the test program**

The movable plunger must be in neutral position. If you are calibrating the unit it is best to turn the stroke adjustment knob so the jaws are open the widest.

Set restart switch to 0.

Press STOP key and hold down while START key is pressed.

Release both keys at the same time.

The message "HT 300 xx.xx" is printed out. (The number xx.xx is important for service, requests etc.).

### **Zero balance**

After starting the test program set the RESTART switch to 1. Press STOP key.

The message "yyyy" is printed out. V



The value "yyyy" must lie between 20 and 60. Is this not the case, adjust the potentiometer "ZERO" so that the value does reach this range. By pressing the STOP-key the zero point can be printed out as often as required. Difference of less than 7 between the single printouts are normal and do not influence the accuracy of the instrument.

### **Calibration**

a) Turn instrument to the side b) Set RESTART-switch on position 2 c) Press STOP-key d) Set RESTART-switch on position 3 e) Place calibration weight

f) Press STOP-key

The message "zzzz N" is printed out (printout always in N, regardless of the setting of the units-switch)

g) Remove calibration weight

h) If necessary adjust potentiometer "GAIN"

1) Repeat operations b) to h) until the balance is satisfactory.

### **Explanation of error information**

NO SAMPLE: No specimen has been inserted.

STOP:

> MAX:                    T.e specimen withstood 300 N.

ZERO: - when pressing the START-key:

The specimen is already pressing on the load cell,  
or: zero point adjustment is incorrect (see test program: zero point adjustment)

- during pre-run:

When pressing the START-key the specimen has pressed stronger on the load cell than during the pre-run.

The STOP-key interrupts the run of the program.

## **CONNECTION POSSIBILITIES**

### **Connection of an external start-key**

It is possible to connect an external start-key (e.g. footswitch) to the 5-pin socket; the function of an external start-key is identical to that of the internal one.

When starting, the sockets 1 and 3 must be short-circuited (30 V, 30 mA).

Socket 5 is earthed when the movable plunger is in its neutral position. (Open-collector-output NPN, allowable load 30 V, 50 mA).

### **Connection to parallel-interface (94-00442) (OPTION)**

.. P

Data from the PTB300 are only given out, if the lead FO is held on OV.

- If the lead ]7S is held on OV, the PTB300 keeps the data output until the lead 13[ also indicates OV.
- The data are valid when the "LOAD"-impulse has started and remain unchanged until the break of the next specimen has occurred.

### **Plug connections**

#### **a) 5-pin socket:**

Pin 1: START-signal (input)

2: no connection

3: earth

4: no connection

5: neutral position of plunger (output)

b) 25-pin socket on parallel interface (TRW-Cinch DBC-25S)

50-pin socket on connecting cable (A~phenol 57-30500)

The program version 8245 of the aTB300 gives following additional possibilities.

- Self test
- printout of result with average value
- equipment optionally with connection for:
  - parallel-interface (PCB 94-00422), qualified for HP 97S
  - serial-interface RS 232-C (PCB 94-00423)

#### SELF TEST

The memory and the electronics are automatically tested when the power is turned on. "ERROR" will be printed when a fault is detected. The unit may still be partly operational. Please inform the service-office.

#### PRINTOUT OF RESULT

With adequate operation the ?TB300 can print out the following information additionally to a measuring series:

- no of samples (batch sizes max. 99 samples or max. 4800 N, 700 SC, 489 kP)
- average value

lowest value

- highest value

#### OPERATION

Selection of the additional printout

The instrument is ready for measurement

- set RESTART-switch to position 9
- press STOP-key and keep it pressed
- press START-key shortly
- release STOP-key

sum of breaking-value batch size

MEASUREMENT is executed normally (with or without RESTART)

### CALLING

the additional information after the last measurement

- set RESTART-switch to position 9
- press STOP-key and keep it pressed
- press START-key shortly
- release STOP-key

With the call of the additional information the selection for the following measuring series is executed at the same time. This means that the measurement for the following series can be done directly, without additional selection.

### DATA-PRINTOUT

NR	17
XM	4.1 KP
XMIN	3.7 KP
XMAX	4.8 KP
17	3.8 KP
16	4.3 KP
15	4.1 KP
14	3.7 KP
13	3.9 KP
12	4.3 KP
11	4.0 KP
^10	4.0 KP
9	4.1 KP
8	4.0 KP
7	3.7 KP

6	4.5 KP
5	3.7 KP
4	3.9 KP
3	4.3 KP
2	4.8 KP
1	4.0 KP

- no. of samples

- average value

- lowest value

- **highest value**

- last measurement

The additional printout ensues automatically if the extreme values are reached in a measuring series:

- after 99 samples

- the sum of the measuring values exceeds 4800 N, 700 SC, 489 kp

INTERFACE RS 232-C

Control switch on the PC/board (switch only when instrument is switched off)

BAUD Rate

STOP BITS

Parity

Parity

Data-Bits

DATA OUTPUT

z

110, 150, 300, 600, 1~00, 2400, 4800, 9600

1

OFF (aus)

ODD (ungerade)

.

7

1

ON (e~n;

EVEN (gerade)

8

**The same information which is output to the printer, is also available on th~ connector.**

DATA INPUT

CT0.N

M

?

1

ACCEPTANCE TIMING

As START-Switch

As STOP-Switch

Continuously

Continuously

Continuously

1 Continuously (except tool is moving fore)

MEAN VALUE

REPEAT OUTPUT

FUNCTION

KEYBOARD

START , START switch

STOP STOP switch

1

Printer ON

—

i Printer OFF

Printer of. PTB300 ON

Printer of PTfi 300 OFF (Except for ERRORmessage)

Sequence of:

RESTART = 9, STOP,START

The previous data output is repeated

Other signs are ignored

If a sign is accepted it is sent back, except the "?", it

directly gives the data output

Is a parity error detected or an overflow of the input

register, then a "?" is sent back.

## **Spare Parts List**

<u>Part Number</u>	<u>Description</u>
15-156-	Power Button
15-156-	Start/Stop Button
15-156-	Button Lamp
15-156-	Button Lens
15-156-	Button Cover Top
15-156-	Button Cover Bottom
15-156-	Restart Switch (Digit Switch)
15-156-	Standard Jaw Set
15-156-0053	Printer Paper