



Key International Inc. 480 Rt. 9 Englishtown, NJ 07726
(732) 536-9700 Phone (732) 972-2630 Fax
www.keyinternational.com

May 30, 2007

HT Calibration Statement

Please be advised that the Key International, Inc. Model HT-300 and HT-500, as well as the HT-300S and HT-500S Tablet Hardness Testers can only be calibrated in the unit of Newtons. This is due to the firmware programming in the unit's EPROM. This unit cannot be changed. When the Calibration function is activated, it automatically sets the unit in Newtons. When the Calibration function is completed, the unit will resume reading in the previously selected units. The proper calibration procedures can be found at the Key International website: www.keyinternational.com.

Signed,

A handwritten signature in black ink, appearing to read "Erik Voorhees". The signature is stylized and cursive.

Erik Voorhees
Electronics Technician
Key International, Inc.
erik@keyinternational.com