Servo R&D Softgel Machine SRD

FEATURES

Changsung is introducing a brand new machine which is servo RND. It has capable to produce both animal and oluscap softgel.





SPECIFICATION

Dimension: 1203(L) × 905(W) × 1622(H) Weight: 670Kg approx Electric power: AC 220, 380, 400V / 50,60 Hz / 3Ph Power consumption: 3.7Kw approx Die roll speed: 3rpm Compressed Air: 5Kg f/m² Die roll Size: Ø100 × 100mm

Special features

Anodizing alumium and stainless steel
Contact part : STS 304 or STS 316
In compliance with GMP
complying with TIR securely
AGT
Stainless steel spreader box
New pinch valve system





1. Die Roll Housing

- Capsule forming device
 - Capsule forming with wedge & die roll combination
- Die roll pressure check system (Load-cell for front and rear)

2. AGT(Automatic Gelatin Thickness)

- This system is interlocked with spreader box
- No need manual thickness gauge tool
- Easy to set up the thickness of gelatin
- Easy device detachable from spreader box for cleaning
- Contol automatically thickness by laser

3. Medichine Pump & Block

- Medichine injection
- Very precise medicine injection by cam
- Titanium coated plungers
- Control medicine weight by touch screen (servo motor driven)

4. Die Roll

Automatic alignment & injection timing by touch screen

5. Spreader Box

- Gelatin ribbon forming device
- From liquid status into solidified riboon
- Automatic gelatin level sensor in spreader box
- AGT conneted

6. Conveyor

Detachable conveyor for easy cleaning

7. Wedge

Injection up& down by manual











1 2 3 4 5 6

7







99, Soheul-Ro, Soheul-Eup, Pocheon-Si, Gyeonggi-Do, Korea TEL 82-31-542-6336 FAX 82-31-542-6543 E-mail master@softgel.kr www.softgel.kr



4 Corporate Drive Cranbury, NJ 08512

P: 609-619-3685 F: 609-619-3686 sales@keyinternational.com www.keyinternational.com